

Product Datasheet

Bovine BMPR1A ELISA Kit (orb438224)

Description	Bovine Bone Morphogenetic Protein Receptor 1A (BMPR1A) Sandwich ELISA Kit
Reactivity	Bovine
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	F1MCE7
Sensitivity	0.53ng/mL
Entrez	404101
Expiration Date	Please enquire.