

Product Datasheet

Human BMPR1B ELISA Kit (orb438225)

Description	Human Bone Morphogenetic Protein Receptor 1B (BMPR1B) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O00238
Sensitivity	0.134ng/mL
Entrez	658
Expiration Date	Please enquire.