

Product Datasheet

Mouse BRAK ELISA Kit (orb438236)

Description	Mouse Breast And Kidney Expressed Chemokine (BRAK) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9WUQ5
Sensitivity	0.054ng/mL
Entrez	57266
Expiration Date	Please enquire.