

Product Datasheet

Human BRCA1 ELISA Kit (orb438238)

Description	Human Breast Cancer Susceptibility Protein 1 (BRCA1) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P38398
Sensitivity	15.6pg/mL
Entrez	672
Expiration Date	Please enquire.