

Product Datasheet

Human BSP ELISA Kit (orb438244)

Description	Human Bone Sialoprotein (BSP) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P21815
Sensitivity	11.2pg/mL
Entrez	3381
Expiration Date	Please enquire.