

Product Datasheet

Mouse bTC ELISA Kit (orb438249)

Description	Mouse Betacellulin (bTC) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q543J8
Sensitivity	7.3pg/mL
Entrez	12223
Expiration Date	Please enquire.