



Product Datasheet

CD40 antibody (orb43825)

Description	Mouse Monoclonal to CD40.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human CD40a
Target	CD40
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The purified antibody is effective for indirect immunofluorescence staining of human cells by flow cytometry analysis using 1 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P05540
Entrez	958
Dilution Range	Flow cytometry: The purified antibody is effective for indirect immunofluorescence staining of human cells by flow cytometry analysis using 1 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human CD40 positive lympho...



Flow cytometry surface staining pattern ...