

Product Datasheet

Mouse C1orf85 ELISA Kit (orb438259)

Description	Mouse Chromosome 1 Open Reading Frame 85 (C1orf85) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9JHJ3
Sensitivity	0.123ng/mL
Entrez	56700
Expiration Date	Please enquire.