

Product Datasheet

Human C1qA ELISA Kit (orb438260)

Description	Human Complement Component 1, Q Subcomponent A (C1qA) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	Assay Length: 1-4.5 hoursMethod: sandwichserum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P02745
Sensitivity	0.059ng/mL
Entrez	712
Expiration Date	Please enquire.