

Product Datasheet

Human C1qL1 ELISA Kit (orb438263)

Description	Human Complement Component 1, Q Subcomponent Like Protein 1 (C1qL1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	075973
Sensitivity	0.053ng/mL
Entrez	10882
Expiration Date	Please enquire.