

Product Datasheet

Rat Complement Component 3 ELISA Kit (orb438280)

Description	Rat Complement Component 3 (C3) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P01026
Sensitivity	10.62ng/mL
Entrez	24232
Expiration Date	Please enquire.