

Product Datasheet

Human C4b ELISA Kit (orb438284)

Description	Human Complement Component 4b (C4b) Sandwich ELISA Kit
Reactivity	Human
Range	4.687-300ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P0C0L5
Sensitivity	1.8ng/mL
Entrez	721
Expiration Date	Please enquire.