

Product Datasheet

Mouse C4b ELISA Kit (orb438285)

Description	Mouse Complement Component 4b (C4b) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01029
Sensitivity	0.245ng/mL
Entrez	12268
Expiration Date	Please enquire.