

Product Datasheet

Human C4BPα ELISA Kit (orb438286)

Description	Human C4 Binding Protein Alpha (C4BPα) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P04003
Sensitivity	0.62ng/mL
Entrez	722
Expiration Date	Please enquire.