

Product Datasheet

Human C4BPb ELISA Kit (orb438287)

Description	Human C4 Binding Protein Beta (C4BPb) Sandwich ELISA Kit
Reactivity	Human
Range	0.468-30ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P20851
Sensitivity	0.17ng/mL
Entrez	725
Expiration Date	Please enquire.