

Product Datasheet

Human Complement Component 9 ELISA Kit (orb438310)

Description	Human Complement Component 9 (C9) Sandwich ELISA Kit
Reactivity	Human
Range	0.468-30ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P02748
Sensitivity	0.19ng/mL
Entrez	735
Expiration Date	Please enquire.