

Product Datasheet

Human CA1 ELISA Kit (orb438317)

Description	Human Carbonic Anhydrase I (CA1) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P00915
Sensitivity	0.58ng/mL
Entrez	759
Expiration Date	Please enquire.