

Product Datasheet

Monkey CA2 ELISA Kit (orb438321)

Description	Monkey Carbonic Anhydrase II (CA2) Sandwich ELISA Kit
Reactivity	Monkey
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6UIR8
Sensitivity	1.22ng/mL
Entrez	574263
Expiration Date	Please enquire.