

Product Datasheet

Human CA4 ELISA Kit (orb438323)

Description	Human Carbonic Anhydrase IV (CA4) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P22748
Sensitivity	0.31ng/mL
Entrez	762
Expiration Date	Please enquire.