

Product Datasheet

Human CA5B ELISA Kit (orb438325)

Description	Human Carbonic Anhydrase VB, Mitochondrial (CA5B) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9Y2D0
Sensitivity	0.058ng/mL
Entrez	11238
Expiration Date	Please enquire.