

Product Datasheet

Human CA6 ELISA Kit (orb438326)

Description	Human Carbonic Anhydrase VI (CA6) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, saliva and other biological fluids.
Sample Types	serum, plasma, saliva and other biological fluids.
Uniprot ID	P23280
Sensitivity	29.3pg/mL
Entrez	765
Expiration Date	Please enquire.