

Product Datasheet

Rat CA9 ELISA Kit (orb438330)

Description	Rat Carbonic Anhydrase IX (CA9) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	F7FL00
Sensitivity	16.2pg/mL
Entrez	313495
Expiration Date	Please enquire.