

Product Datasheet

Human CACNa1B ELISA Kit (orb438331)

Description	Human Calcium Channel, Voltage Dependent, N-Type, Alpha 1B Subunit (CACNa1B) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q00975
Sensitivity	0.129ng/mL
Entrez	774
Expiration Date	Please enquire.