

Product Datasheet

Human CACNa1D ELISA Kit (orb438334)

Description	Human Calcium Channel, Voltage Dependent, L-Type, Alpha 1D Subunit (CACNa1D) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q01668
Sensitivity	0.121ng/mL
Entrez	776
Expiration Date	Please enquire.