

Product Datasheet

Rat CACT ELISA Kit (orb438340)

Description	Rat Carnitine Acylcarnitine Translocase (CACT) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P97521
Sensitivity	0.087ng/mL
Entrez	117035
Expiration Date	Please enquire.