

Product Datasheet

Human CALB ELISA Kit (orb438345)

Description	Human Calbindin (CALB) Sandwich ELISA Kit
Reactivity	Human
Range	12.5-800pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P05937
Sensitivity	5.7pg/mL
Entrez	793
Expiration Date	Please enquire.