

Product Datasheet

Human CAMK2b ELISA Kit (orb438357)

Description	Human Calcium/Calmodulin Dependent Protein Kinase II Beta (CAMK2b) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q13554
Sensitivity	0.113ng/mL
Entrez	816
Expiration Date	Please enquire.