

Product Datasheet

Rat CAMK2g ELISA Kit (orb438358)

Description	Rat Calcium/Calmodulin Dependent Protein Kinase II Gamma (CAMK2g) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P11730
Sensitivity	0.123ng/mL
Entrez	171140
Expiration Date	Please enquire.