

Product Datasheet

Human CAMKK2 ELISA Kit (orb438359)

Description	Human Calcium/Calmodulin Dependent Protein Kinase Kinase 2 (CAMKK2) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q96RR4
Sensitivity	0.114ng/mL
Entrez	10645
Expiration Date	Please enquire.