

Product Datasheet

Mouse CAPN10 ELISA Kit (orb438366)

Description	Mouse Calpain 10 (CAPN10) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9ESK3
Sensitivity	0.103ng/mL
Entrez	23830
Expiration Date	Please enquire.