

Product Datasheet

Mouse CAPN1 ELISA Kit (orb438368)

Description	Mouse Calpain 1, Large Subunit (CAPN1) Sandwich ELISA Kit
Reactivity	Mouse
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	O35350
Sensitivity	0.58ng/mL
Entrez	12333
Expiration Date	Please enquire.