

Product Datasheet

Rat CART ELISA Kit [Out of stock] (orb438378)

Description	Rat Cocaine And Amphetamine Regulated Transcript (CART) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P49192
Sensitivity	5.8pg/mL
Entrez	29131
Expiration Date	Please enquire.