

Product Datasheet

Rat CASP14 ELISA Kit (orb438381)

Description	Rat Caspase 14 (CASP14) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	D3ZZ65
Sensitivity	0.103ng/mL
Entrez	299587
Expiration Date	Please enquire.