

Product Datasheet

Mouse CAST ELISA Kit (orb438403)

Description	Mouse Calpastatin (CAST) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P51125
Sensitivity	0.043ng/mL
Entrez	12380
Expiration Date	Please enquire.