

Product Datasheet

Rat CAT2 ELISA Kit (orb438405)

Description	Rat Cationic Amino Acid Transporter 2 (CAT2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	B5D5N9
Sensitivity	0.069ng/mL
Entrez	64554
Expiration Date	Please enquire.