

Product Datasheet

Chicken CAV1 ELISA Kit (orb438410)

Description	Chicken Caveolin 1 (CAV1) Sandwich ELISA Kit
Reactivity	Gallus
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P35431
Sensitivity	0.064ng/mL
Entrez	373996
Expiration Date	Please enquire.