

Product Datasheet

Goat CAV1 ELISA Kit (orb438411)

Description	Goat Caveolin 1 (CAV1) Sandwich ELISA Kit
Reactivity	Goat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q069J5
Sensitivity	0.06ng/mL
Entrez	100860749
Expiration Date	Please enquire.