

Product Datasheet

Mouse CAV1 ELISA Kit (orb438413)

Description	Mouse Caveolin 1 (CAV1) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P49817
Sensitivity	10.7pg/mL
Entrez	12389
Expiration Date	Please enquire.