

Product Datasheet

Mouse CbS ELISA Kit (orb438422)

Description	Mouse Cystathionine Beta Synthase (CbS) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q91WT9
Sensitivity	0.059ng/mL
Entrez	12411
Expiration Date	Please enquire.