

Product Datasheet

Mouse CBX3 ELISA Kit (orb438424)

Description	Mouse Chromobox Homolog 3 (CBX3) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P23198
Sensitivity	0.046ng/mL
Entrez	12417
Expiration Date	Please enquire.