

Product Datasheet

Human CC16 ELISA Kit (orb438425)

Description	Human Clara Cell Protein 16 (CC16) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, lung lavage fluid , cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, lung lavage fluid , cell culture supernatants and other biological fluids.
Uniprot ID	P11684
Sensitivity	23.5pg/mL
Entrez	7356
Expiration Date	Please enquire.