

Product Datasheet

Rat CC16 ELISA Kit (orb438427)

Description	Rat Clara Cell Protein 16 (CC16) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, lung lavage fluid , cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, lung lavage fluid , cell culture supernatants and other biological fluids.
Uniprot ID	P17559
Sensitivity	35pg/mL
Entrez	25575
Expiration Date	Please enquire.