

Product Datasheet

Human CC2D1A ELISA Kit (orb438428)

Description	Human Coiled Coil And C2 Domain Containing Protein 1A (CC2D1A) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q6P1N0
Sensitivity	0.128ng/mL
Entrez	54862
Expiration Date	Please enquire.