

Product Datasheet

Human CCbL1 ELISA Kit (orb438429)

Description	Human Cysteine Conjugate Beta Lyase, Cytoplasmic (CCbL1) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q16773
Sensitivity	27pg/mL
Entrez	883
Expiration Date	Please enquire.