



Product Datasheet

CD42b antibody (orb43843)

Description	Mouse Monoclonal to CD42b.
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr
Immunogen	Peripheral blood mononuclear cells of a patient suffering with CLL.
Target	CD42b
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P07359
Entrez	2811
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human thrombocytes (red-fi...



Flow cytometry surface staining pattern ...