

Product Datasheet

Human CCDC80 ELISA Kit (orb438431)

Description	Human Coiled Coil Domain Containing Protein 80 (CCDC80) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q76M96
Sensitivity	0.053ng/mL
Entrez	151887
Expiration Date	Please enquire.