

Product Datasheet

Canine CCKAR ELISA Kit (orb438437)

Description	Canine Cholecystokinin A Receptor (CCKAR) Sandwich ELISA Kit
Reactivity	Canine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q5D0K2
Sensitivity	0.124ng/mL
Entrez	488846
Expiration Date	Please enquire.