

Product Datasheet

Human CCKAR ELISA Kit (orb438440)

Description	Human Cholecystokinin A Receptor (CCKAR) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P32238
Sensitivity	0.084ng/mL
Entrez	886
Expiration Date	Please enquire.