

Product Datasheet

Rabbit CCKAR ELISA Kit (orb438443)

Description	Rabbit Cholecystokinin A Receptor (CCKAR) Sandwich ELISA Kit
Reactivity	Rabbit
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O97772
Sensitivity	0.61ng/mL
Entrez	100009241
Expiration Date	Please enquire.