

Product Datasheet

Rat CCK ELISA Kit (orb438449)

Description	Rat Cholecystokinin (CCK) competitive inhibition ELISA Kit
Reactivity	Rat
Range	12.5-800pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01355
Sensitivity	3.75pg/mL
Entrez	25298
Expiration Date	Please enquire.